

Dissection of human and mouse embryonic lungs and set-up of organoid culture

ER Emma L Rawlins

Updated date: Sep 28, 2021

 An abbreviated version of this protocol was published in eLIFE in Jun 2017

Human embryonic lung epithelial tips are multipotent progenitors that can be expanded in vitro as long-term self-renewing organoids

DOI: 10.7554/eLife.26575

Detailed protocol

Human fetal lung dissection

- Collect lungs in Hibernate A medium & keep it at 4°C
- On the dissection day need to:
 - Thaw Matrigel in the back of the tissue culture room fridge (takes 3 hr to thaw)
 - Place a 48-well plate (greiner bio-one, Cellstar, Cat No 677102) in the big incubator to warm up the plates, so that the Matrigel will set quickly and form a mound instead of spreading all over the well.
- Wash tissues with T/C PBS.
- Use forceps to remove esophagus and trachea.
- Treat lung with dispase (8 U/ml PBS) at room temperature for 2 min. (Cut lungs into smaller pieces if big.)
- Wash off dispase with PBS & keep the lungs in Advanced+++ DMEM/F12 on ice.
- Dissect a small piece of lung under microscope in a 5 cm petri-dish containing some Advanced+++ DMEM/F12; remove mesenchyme to reveal tips with tungsten needles.
- Cut tips with a sharp tungsten needle on the bud necks.
- Transfer tips with a 50 ul pipette to another petri-dish containing some Advanced+++ DMEM/F12.
- At the same time, prepare cold Matrigel in Eppendorf tubes in a T/C hood.
 - Cool Matrigel and 1.5ml Eppendorf tubes in an ice box.
 - 30ul Matrigel per Eppendorf tube.
- Group 5 tips together and pipette them to a Matrigel Eppendorf prepared above with a 50 ul pipette under microscope. Keep on ice.
- After all the tips have been transferred to the Matrigel eppendorfs, pipette the tips containing Matrigel above to a 48-well plate.
 - Avoid bubbles
 - ~ 25 ul each well
 - ~ 5 tips per well
- Place the 48-well plate in a 37°C incubator for 5 min to solidify the Matrigel
- Add self renewal medium along the side of the well (300 ul per well).
- Fill empty wells with PBS to minimise evaporation from Matrigel wells.
- Incubate in a 37°C CO2 incubator.
- Change medium on Tuesday and Friday; prepare fresh medium weekly
- Passage:
 - 1:4 dilution for the first passage, and then 1:5 dilution for further passages

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Rawlins, E. (2021). Dissection of human and mouse embryonic lungs and set-up of organoid culture. Bio-protocol Preprint. bio-protocol.org/prep1383.
2. Nikolić, M. Z., Caritg, O., Jeng, Q., Johnson, J., Sun, D., Howell, K. J., Brady, J. L., Laresgoiti, U., Allen, G., Butler, R., Zilbauer, M., Giangreco, A. and Rawlins, E. L. (2017). Human embryonic lung epithelial tips are multipotent progenitors that can be expanded in vitro as long-term self-renewing organoids. eLIFE. DOI: [10.7554/eLife.26575](https://doi.org/10.7554/eLife.26575)

Copyright: Content may be subjected to copyright.

